



P.O. Box 746, 5340 AS Oss +  
The Netherlands +  
T +31 412 67 13 00 +  
F +31 412 67 13 10 +  
[www.medica-europe.nl](http://www.medica-europe.nl) +



## THERE ARE TWO TYPES OF MEDICA BRUSH®

### DRY STERILE (art. nr. 255A480)

- sterilized using ethyleneoxide (ETO)
- without soap
- easy to open
- can also be used for cleaning the intact patients skin

### IMPREGNATED WITH CHLORHEXIDINEDIGLUCONATE IN DETERGENT (art. nr. 255A481)

- contains 4% chlorhexidinedigluconate
- broad spectrum of action
- good residual working
- mild cleansing agent: minimal chance of skin irritations
- enough soap for 3-5 minutes' washing
- easy to open

**MEDICA**BRUSH®

**INFECTION PREVENTION  
WITH A SOFT TOUCH**



INFECTION PREVENTION  
WITH A SOFT TOUCH

A good hygiene forms the basis of preventing infection in the hospital. With good reason, since the majority of hospital infections are the result of insufficient hand and skin hygiene. The introduction of the Medica Brush® in your hospital will enable you to establish the right conditions for optimal preoperative hygiene of hand and skin.

The disposable Medica Brush® has many advantages which traditional reusable brushes simply cannot offer. The Medica Brush® is very convenient to use. Its soft bristles for nails, knuckles and skin clean effectively while being kind to the skin. The attached sponge gently cleans the hands, forearms and skin avoiding the risks and discomfort of skin irritation and abrasion. A separate nail cleaner is also enclosed.

The fact that brush, soap and nail cleaner are combined in one package means that the Medica Brush® also reduces handling costs. All these features combined make Medica Brush® the complete solution for good preoperative hygiene of hand and skin.

The Medica Brush® comes in a wellsealed package which keeps the soap stable but remains very simple and hygienic to open. The handy dispenser box contains 25 brushes and can be attached to the wall beside the wash basin and/or in the OR. A holder is available in a choice of aluminium or transparent plastic. The strictest of quality controls (ISO 13485) are applied throughout the production process and the Medica Brush® complies fully with the very latest standards.

PHYSICAL AND CHEMICAL DETAILS

	Chlorhexidine
pH value	+/- 6,0
Density (g/cm³)	+/- 1,02
Appearance	amber liquid

All the benefits of the disposable Medica Brush® at a glance:

- significant reduction of micro-organisms (Medica Brush® chlorhexidine)
- deeper-lying resident bacterial flora is not removed unnecessarily
- soft and kind to the skin, encouraging more thorough use
- effective, sensitive and conveniently shaped brush, comprising two layers of ultrasoft material
- balanced supply of soap through the squeeze-through baseplate
- no waste of soap
- no need for resterilization procedures, saving time, work and money



THE MEDICA BRUSH® WORKS FOR  
HANDS AS FOLLOWS:



Take a Medica Brush® from the dispenser and open the packaging.



Wet your hands and forearms under the tap. Use the sponge to wash your hands and forearms.



Clean your nails with the accompanying nail cleaner.



Use the bristles of the Medica Brush® to clean your nails and knuckles.



Carefully rinse and dry your hands and forearms.

Recommendation: After using the Medica Brush®, apply a hand disinfection.